

Natural Rubber Latex Powder Free Gloves – Extended Cuff High Risk, Blue

User Instructions

1.0 Description of the device:

Natural Rubber Latex Powder Free Gloves – Extended Cuff High Risk, Blue meet all the requirements of ASTM D3578-19 Standard Specification for Rubber Examination Gloves.

2.0 Indication for Use of the Device

Natural Rubber Latex Powder Free Gloves – Extended Cuff High Risk, Blue which is a disposable device intended for medical purpose that is worn on the examiner's hand or finger to prevent contamination between examiner and patient.

3.0 Labelling

The user instructions are demonstrated on the actual packaging material which will be at the final packaging of goods.

4.0 Pre-Use Information

4.1 Safety Considerations

If reaction in the form of skin irritation is noticed, remove gloves immediately and wash affected part with water. If there is no relief, seek medical reactions.

4.2 Limitations of use

Natural Rubber Latex Powder Free Gloves – Extended Cuff High Risk, Blue is intended for one use, or for use a single patient during a single procedure.

4.3 Marking Recommendation and Restriction

- 4.3.1 Resistance against Bacteria and Fungi- Pass
- 4.3.2 Resistance against Virus- Pass
- 4.3.3 This glove is for single use only!
- 4.3.4 This emergency medical protective single-use examination glove meets requirements of NFPA 1950 (1999), 2025 edition.
- 4.3.5 Do not remove this label!
- 4.3.6 Protective single use examination gloves for emergency medical operations in accordance with NFPA 1950-2025.
- 4.3.7 The penetration resistance has been assessed under laboratory condition and relates only to related specimens.
- 4.3.8 Before usage, inspect gloves for any defect or imperfections. Please make sure gloves are free from defects or abnormality such as tear, hole and pinhole. If found any gloves with defects do not use the gloves. Get new gloves. Ensure chemicals cannot enter via cuff. Remove the glove immediately if contaminated by a concentrate spill.
- 4.3.9 It is recommended to check that the gloves are suitable for the intended use because the conditions at the workplace may differ from the type of test depending on temperature, abrasion and degradations.
- 4.3.10 When used, protective gloves may provide less resistance to the dangerous chemical due to changes in physical properties. Movements, snagging, rubbing, degradation caused by the chemical contact ect. may reduce the actual use time significantly.

4.3.11 Warning

This product is not made for fire protection use. This product is non-sterile and single use only.

4.3.12 Caution

This product contains Natural Rubber Latex which may cause allergic reactions including anaphylactic responses in some individuals, and latex sensitive individuals. Discontinue use immediately and consult physician if allergic reactions occur.

4.4 Statement

4.4.1 “Most performance properties cannot be tested by the user in the field.”

4.5 Warranty information

Natural Rubber Latex Powder-Free Gloves – Extended Cuff, High Risk, Blue is a single-use device. Recipients of our products are responsible for complying with all applicable laws and regulations. This warranty is valid for five years from the date of manufacture, provided that the gloves have been stored in their original packaging, unopened packaging, kept in a cool, dry, and clean environment, and not exposed to UV light, excessive heat, or humidity. This warranty applies only to unused and unmarked gloves. Once the gloves are removed from their packaging or used, the warranty becomes void.

5.0 Preparation for Use

5.1 Sizing/adjustment

Sizing measurements are demonstrated on the actual packaging material which will be at the final packaging of goods.

5.2 Recommended storage practices

Store in a dry, cool and ventilated area. Shield open box from direct sunlight, fluorescent lighting and x-rays. Improper storage will decrease usable life.

6.0 Inspection frequency and details

6.1 Before usage, inspect gloves for any defect or imperfections. Please make sure gloves are free from defects or abnormality such as tear, hole and pinhole. If found any gloves with defects do not use the gloves. Get new gloves. Ensure chemicals cannot enter via cuff. Remove the glove immediately if contaminated by a concentrate spill.

7.0 Don/doff

7.1 Donning and doffing procedures

7.1.1 Donning Instructions

Inspect the glove and ensure there is no pinholes and tears. Insert all five fingers into the cuff and pull the cuff over your wrist until its properly place. Check that the glove's fit is secure around the fingers and the palm. Also check the cuff, which should have a snug fit around your wrist. If the fit feels too tight or too loose, consider changing size to avoid any tears or discomfort.

7.1.2 Doffing Instructions

Using a gloved hand, Grasp the outside edge of the glove near the wrist, be careful to not touch your skin. Peel the glove away from the hand,

turning it inside out. Hold it in the opposite gloved hand. Slide an unglove finger under the wrist of the remaining glove, being careful not to touch the outside of the glove. Peel the remaining glove off from the inside, creating a “bag” containing both gloves.

7.2 Sizing and adjustment procedure

- 7.2.1 Gloves are available in multiple sizes and must be properly fitted to ensure protection, dexterity, and full coverage.

7.3 Interface issues

- 7.3.1 When used, protective gloves may provide less resistance to the dangerous chemical due to changes in physical properties. Movements, snagging, rubbing, degradation caused by the chemical contact etc. may reduce the actual use time significantly.

8.0 Proper Use consistent with NFPA 1500, NFPA 1581, 29 CFR 1910.132, “General Requirements of subpart I, Personal Protective Equipment,” and 29 CFR 1910.1030, “Protecting Health Care Workers from Occupational Exposure to Blood-Borne Pathogens”

9.0 Decontamination procedure

- 9.1 This product is non-sterile and single use only. Removed carefully to avoid contact with the outer surface. Do not wash, disinfect, or reuse. Hand hygiene must be performed after glove removal using soap and water or an appropriate alcohol-based hand rub.

10.0 Disposal Criteria and Considerations

10.1 Disposal

Dispose of in accordance with all local, state and federal regulations. All empty packaging should be disposed of in accordance with local, state, and federal facility regulations or recycled/reconditioned at an approved facility.